## ORBIS Acryl Xject

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### EN Processing instructions

#### 1 Intended purpose

Denture base materials consist of a liquid and a powder component. A dental technician uses these materials to produce splints or dentures in different production processes at the dental lab. Denture base materials are intended for the fabrication of splints, partial dentures, full dentures and model dentures.

### 2 Description of product and users

#### 2.1 Product description

ORBIS Acryl Xject is a cold-curing methylmethacrylate-based denture acrylic for the fabrication of dentures using various manufacturing processes.

### 2.2 Users

For use by laboratory technicians in a dental laboratory

### 3 Composition

- Powder:** Pearl polymerisate made from polymethyl methacrylate (PMMA) pigments, initiators
- Liquid:** Methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilisers
- ORBIS Acryl Xject is cadmium-free.

### 4 Indications

- For the production of partial and full dentures using the injection, flow or flask technique
- Production of the acrylic portion of implant-borne dental appliances and model casting dentures
- Repairs and relining of existing restorations.
- For the production of splints
- Mouthguards
- Bite splints

### 5 Contraindications

- Do not use in the case of a known allergy to one of the components.

### 6 Warnings

- Liquid:** Contains methyl methacrylate, tetramethylene dimethacrylate.
- Danger. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation.
- Powder:** Contains dibenzoyl peroxide, methyl methacrylate. May produce an allergic reaction. Harmful to aquatic life with long lasting effects.

### 7 Safety instructions

- Liquid:** Keep away from heat/sparks/open flames/hot surfaces. No smoking Use explosion-proof electrical/ventilating/lighting/equipment. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/shower.
- Powder:** Avoid release to the environment.

### 8 Interactions with other medicinal products

The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

### 9 Application / Preparation

#### 9.1 Mixing ratio

- 2 units powder to 1 unit liquid by weight
- On average, partial dentures require: 30 g powder: 15 g liquid (weight)
- Veins can be included, if desired. Please observe separate manufacturer's information and instructions.

### 9.2 Production of partial and full dentures as well as implant-borne appliances and model casting completions

#### 9.2.1 Preparatory work

Produce the plaster cast and wax-up of the teeth. The fabricated tooth set in wax is secured with plaster or silicone (matrix/flask). Then the wax is eliminated. To prevent artificial teeth from breaking off, they must be roughened, have undercuts made, and be wetted with monomer or a suitable bonder. The surfaces of the plaster models are separated with an alginate separating agent according to the manufacturer's instructions.

### 9.2.2 Production using the flow and flask techniques

Mix powder and liquid (according to mixing ratio) and fill into the mould as usual. After mixing, ORBIS Acryl Xject is pourable for approx. 3 minutes at a room temperature of 21 °C; after 4-5 minutes, it transitions into the plastic phase and after no more than 7-8 minutes, the dental appliance must be placed in the pressure vessel. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45 °C.

#### 9.3 Production using the UNIPRESS injection technique (observe separate instructions) (Please follow manufacturer's instructions in the case of other suitable injection systems)

Mix powder and liquid (see mixing ratio above) and after allowing them to thicken for 15 sec., homogeneously mix in a porcelain beaker using a spatula and fill into the UNIPRESS cylinder. As soon as a viscoelastic consistency has been reached, the acrylic is injected into the prepared UNIPRESS flask. These work steps must be completed within a maximum of 8 minutes. Polymerisation takes place under the defined pressurisation of the plate spring system in the UNIPRESS. After 20-30 minutes, the ORBIS Acryl Xject dentures can be deflasked.

### 9.3.1 Repairs and relinings of existing acrylic restorations

Create a model or a matrix, depending on the construction. Clean and roughen the acrylic surfaces to be repaired, smooth the edges and wet with monomer or a suitable bonder. Mix the acrylic in accordance with the mixing ratio and pour in. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45 °C.

### 9.3.2 Finishing the material

After polymerisation lift the workpiece off the model and finish with suitable tools. Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide burrs and diamond instruments are suitable for finishing. Sandpaper or pumice stone is used for pre-polishing.

### 9.3.3 Polishing

The material is polished with goat bristle brushes, polishing paste and soft cotton wool discs. Careful surface finishing and polishing is necessary for an optimal result.

### 9.4 Production of splints

#### 9.4.1 Preparatory work

Produce a plaster model (cast) and block out undercuts with suitable materials. Model the splint in wax as usual. Next, make an impression, create flow channels, and remove wax. The surfaces of the plaster models are separated with an alginate separating agent according to the manufacturer's instructions.

### 9.4.2 Processing

Mix powder and liquid (according to mixing ratio) and fill into the mould as usual. These work steps must be completed after no more than 7-8 minutes.

The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45 °C.

### 9.4.3 Finishing the material

After polymerisation lift the workpiece off the model and finish with suitable tools. Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide burrs and diamond instruments are suitable for finishing. Sandpaper or pumice stone is used for pre-polishing.

### 9.4.4 Polishing

The material is polished with goat bristle brushes, polishing paste and soft cotton wool discs. Careful surface finishing and polishing is necessary for an optimal result.

### 9.5 Repairs

Create a model or a matrix, if necessary. Clean and roughen the acrylic surfaces to be repaired, smooth the edges and wet with monomer or a suitable bonder. Mix the acrylic in accordance with the mixing ratio and pour in. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45 °C.

### 10 Notes on processing

To further reduce the risk of intolerance reactions, we recommend placing the denture in lukewarm water for at least 12 hours prior to insertion.

### 11 Troubleshooting / FAQ Liste

Problem	Cause	Remedy
Acrylic does not cure properly	Incorrect mixing ratio	- Refer to the instructions for use <p>- Ensure correct mixing ratio</p>
Acrylic is whitish at base	- Inadequate insulation <p>- Cast not steeped in water</p>	- Ensure adequate insulation <p>- Check insulating materials if necessary</p> <p>- Allow cast to steep for approx. 10 min in lukewarm water</p>
Acrylic forms bubbles at base	Cast not steeped in water	Allow cast to steep for approx. 10 min in lukewarm water
Acrylic surface white/porous	Acrylic injected too late or placed in pressure vessel too late	Observe the processing times
Differences in colouration following repairs	Incorrect mixing ratio	Ensure correct mixing ratio
Acrylic breaks when lifted out	Inadequate insulation	Check insulation
Acrylic/denture has whitish streaks	- Inadequate mixing <p>- Too much powder used</p>	- Carefully stir acrylic dough with spatula <p>- Ensure correct mixing ratio</p>
Acrylic sets too fast	Excessively high temperatures/excessively long work steps	Observe the temperatures/processing times
Splint breaks	Contouring too thin/unsuitable	Ensure sufficient layer thickness (at least 2-3 mm)

### 12 Storage and handling information

Storage temperature 10-25 °C / 50-77 °F. Close container carefully.

### 13 Shelf life

The maximum shelf life is printed on the label of each pack. Do not use after the expiry date.

### 14 Warnings on side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

### 15 Instructions for disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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