**Product description:**  
ORBIS Compomer is a universal light-curing radiopaque compomer restorative material for filling cavities in both the anterior and posterior areas. It allows a simplified and rapid filling technique as well as combines the fluoride release of glass ionomers with the physical durability and high aesthetics of composites. ORBIS Compomer is used in combination with a universal self-etching bond which provides uniformly high adhesion to enamel and dentine.

**Indications:**  

- filling deciduous teeth
- class V fillings (cervical caries, root erosions, V-shaped defects)
- facetting of discoloured anteriors
- class III fillings
- temporary class I and II fillings
- extended fissure sealing

**Contraindications:**  

- If there known allergies to its ingredients (methacrylates, amines)
- If a dry work field cannot be secured.
- If the prescribed application technique cannot be employed.
- For direct and indirect pulp capping
- For load-bearing permanent fillings in secondary dentition.

**Application:**  
**Preparation/shade selection:**  
Clean the teeth to be treated with fluoride-free cleaning paste. To select the shade, compare the shade guide with original material (corresponds to the shade of an applied and thoroughly polymerised filling) with the cleaned tooth before anaesthesia (moisten shade guide and the tooth, select the shade in daylight). Secure a dry working field.

**Cavity preparation:**  
In general, preparation should be minimal (adhesive filling technique) to conserve healthy tooth substance.  
Observe special preparation forms for extended fissure sealing and on deciduous teeth.  
Caries-free cervical defects do not require preparation. Clean the defect with either cleaning paste or pumice using a rubber cup or rotating brush.

**Cleaning:**  

- Remove all residues from the cavity with the water jet.
- Avoid any contamination with blood or saliva after the cleaning.
- Remove excess water with a faint air jet, yet do not overdry the dentine. A moist dentine surface is necessary for wet-bonding.

**Pulp protection:**  
A calcium hydroxide lining should be applied in close proximity to the pulp.

**Application of Bonding:**  
See corresponding instructions for use.

**Application of ORBIS Compomer:**  
Apply and light-cure fillings more than 2 mm deep in layers. Polymerise each layer for 40 s with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm<sup>2</sup> on halogen polymerisation devices and 300 mW/cm<sup>2</sup> on LED devices.

Use transparent strips and crowns or light wedges. Ensure a tight marginal seal. Hold the tip of the light rod as close as possible to the surface of the filling. If the distance exceeds 5 mm the curing depth may be reduced to polymerisation in the range of the light cone only.

Incomplete polymerisation may lead to discoloration and pulpitis-like complaints.

Finishing and polishing of the filling can be conducted with cooling provided (e. g. fine or extra fine finishing diamonds, polishing discs) immediately after removing the shaping aids. As the last step, the margin of the filling or the entire tooth should be fluoridated.

**Warnings, precautionary measures:**  

- ORBIS Compomer contains Bis-GMA, diurethandimethacrylate, TEGDMA, BHT. ORBIS Compomer should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Preparations containing eugenol or oil of cloves disturb the polymerisation of resin-based filling materials and should not be used in combination with resin-based filling materials.
- Contact with cationic mouth rinses as well as plaque indicators and chlorhexidine may lead to discoloration.

**Storage:**  
Store at temperatures of 4 °C to 23 °C. Close syringes and bond bottles immediately after dispensing to avoid premature curing and/or evaporation of the solvent. Do not use after expiry date.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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