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EN **Instructions for use**

1. Intended purpose

Denture base materials are intended for the fabrication of splints, partial dentures, full dentures and model dentures.

2. Description of product and users

2.1 Description of product

Acryl Base Cold is a cold-curing methylmethacrylate-based denture acrylic used in various processes to fabricate dentures. It is available in various colours, allowing perfect mimicking of form and colour.

2.2 Users

Denture base materials are used by laboratory technicians in a dental laboratory.

3. Composition

Powder: Pearl polymerisate made from polymethyl methacrylate (PMMA) pigments, initiators.

Liquid: Methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilisers.

4. Indications

For the fabrication of partial and full dentures using the flow technique.

Repairs, extensions and relinings of existing restorations using the flow technique.

For the production of splints:

- Mouthguards
- Bite splints

5. Contraindications

Do not use in the case of a known allergy to one of the components.

6. Warnings

Liquid: Contains methyl methacrylate, tetramethylene dimethacrylate. Danger. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation.

Powder: Contains benzoyl peroxide. May produce an allergic reaction. Harmful to aquatic life with long lasting effects.

7. Safety instructions

Liquid: Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use explosion-proof [electrical/ventilating/lighting] equipment. Avoid breathing vapours/spray. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].

Powder: Avoid release to the environment.

8. Interactions with other medicinal products

The dentist should consider known interactions and cross-reactions of the medical product with other materials already in the patient's mouth before using the product.

9. Application / Preparation

9.1 Mixing ratio

10 parts by weight powder, 7 parts by weight liquid

Veins can be included, if desired.

Please observe separate manufacturer's instructions.

9.2 Production of partial and full dentures as well as implant-borne appliances and model casting completions using the flow techniques

9.2.1 Preparatory work

Produce the plaster cast and wax-up of the teeth as usual.

The fabricated tooth set in wax is secured with plaster or silicone (matrix/flask). The wax is then eliminated as usual.

To prevent artificial teeth from breaking off, they must be roughened, have undercuts made, and be wetted with monomer or a suitable bonder.

The surfaces of the plaster models are separated with an alginate separating agent according to the manufacturer's instructions.

9.2.2 Production using the flow technique

Mix powder and liquid (according to the mixing ratio) and after allowing it to thicken for approx. 15 seconds, stir with a spatula until a homogeneous mixture is formed. Fill acrylic dough into the mould as usual. After mixing, Acryl Base Cold is pourable for approx. 2 minutes at a room temperature of 21°C; after 3 minutes, it transitions into the plastic phase. After the start of the plastic phase, the acrylic no longer flows out of the matrix and it can be modelled. The dental appliance must be placed in the pressure vessel after no more than 5 minutes.

The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45°C.

If Acryl Base Cold LF liquid is used, the acrylic is pourable for approx. 3 minutes; after 4 minutes, it transitions

into the plastic phase, and after no more than 6-7 minutes, the dental appliance must be placed in the pressure vessel. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2-5 bar and a water temperature of 45°C.

9.3 Production of splints

9.3.1 Preparatory work

Produce a plaster model (cast) and block out with suitable materials. Model the splint in wax as desired. Next, make an impression, create flow channels, and remove wax. Separate the surface of the plaster model using suitable materials according to manufacturer's instructions.

9.3.2 Processing

Mix powder and liquid (according to the mixing ratio) and after allowing it to thicken for approx. 15 seconds, stir with a spatula until a homogeneous mixture is formed. Fill acrylic into the mould as usual. These work steps must be completed within a maximum of 6 minutes. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2-5 bar and a water temperature of 45°C.

9.4 Finishing the material

After polymerisation, the splint is lifted from the model and finished using suitable instruments. Silicone polishes (square-edge or knife-edge wheels or cylinders), tungsten carbide burrs and diamond instruments are suitable for finishing. Pre-polish with sandpaper or pumice powder. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2-5 bar and a water temperature of 45°C.

9.5 Polishing

The material is polished with goat bristle brushes, polishing paste and soft cotton wool discs. Careful surface finishing and polishing is necessary for an optimal result.

9.6 Repairs

Create a model or a matrix, if necessary. Clean and roughen the acrylic surfaces to be repaired, smooth the edges and wet with monomer or a suitable bond. Mix the acrylic in accordance with the mixing ratio and pour in. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2-5 bar and a water temperature of 45°C.

10. Notes on processing

In the event of deviations from 21°C room temperature:

Higher temperatures shorten, and lower temperatures prolong the above times. To further reduce the risk of intolerance reactions, we recommend placing the denture in lukewarm water for at least 12 hours prior to insertion.

11. Troubleshooting / FAQ List

Problem	Cause	Remedy
Acrylic does not cure properly	Incorrect mixing ratio	- Refer to the instructions for use - Ensure correct mixing ratio
Acrylic is whitish at base	- Inadequate insulation - Cast not steeped in water	- Ensure adequate insulation - Check insulating materials if necessary - Allow cast to steep for approx. 10 min in lukewarm water
Acrylic forms bubbles at base	- Cast not steeped in water	- Allow cast to steep for approx. 10 min in lukewarm water
Acrylic surface white/porous	Acrylic injected too late or placed in pressure vessel too late	Observe the processing times
Differences in colouration following repairs	Incorrect mixing ratio	Ensure correct mixing ratio
Acrylic breaks when lifted out	Inadequate insulation	Check insulation
Acrylic/denture has whitish streaks	- Inadequate mixing - Too much powder used	- Carefully stir acrylic dough with spatula - Ensure correct mixing ratio
Acrylic sets too fast	Excessively high temperatures/ excessively long work steps	Observe the temperatures/ processing times
Splint breaks	Contouring too thin/unsuitable	Ensure sufficient layer thickness (at least 2-3 mm)

12. Information on storage and handling

Storage temperature 10-25°C / 50-77°F. Close container carefully.

13. Shelf life

The maximum shelf life is printed on the label of each pack. Do not use after the expiry date.

14. Warnings on side effects

With proper preparation and use of this medical product, adverse effects are extremely rare. Immune reactions (such as allergies) or local discomfort, however, cannot be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

15. Instructions for disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

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MANI MEDICAL GERMANY GmbH