## EN INSTRUCTIONS FOR USE MD EU Medical device

**Product description:**

ORBI-Bond Universal is a dental adhesive for use with methacrylate-based restorative, core build-up or luting materials. It can be applied with phosphoric acid etching (selective enamel etching or total-etch) of the dental hard tissue as well as without (self-etch).

**Indications:**

[A1] Adhesive for direct restorations of all classes of cavities using light-curing composite / compomer / Ormocer® materials with a methacrylate base.

[A2] Adhesive for core build-ups made of light-curing composite materials.

[A3] Adhesive for indirect restorations with light-curing luting composites.

[B] Intraoral repair of composite, compomer or Ormocer® restorations, ceramic veneers and all-ceramic restorations without an additional primer.

[C] Treatment of hypersensitive tooth necks.

[D] Protective varnish for glass ionomer cement restorations.

[E1] Sealing of cavities prior to amalgam restorations.

[E2] Sealing of cavities and core preparations prior to the temporary luting of indirect restorations.

In combination with **ORBI-Bond Universal DC Activator (Dual Cure Activator)**:

[F1] Adhesive for direct self-curing or dual-curing composite restorations of all classes of cavities and for core build-ups.

[F2] Adhesive for indirect restorations using dual-curing or self-curing luting composites.

[G] Luting of root posts with dual-curing or self-curing luting composites.

**Method of use:**

ORBI-Bond Universal not only ensures a secure bond to natural dental hard tissue, but also to composites, metals and ceramics.

It is necessary to isolate the operating site from moisture. Any contamination of the cavity or material surface with blood or saliva must be prevented. The use of a rubber dam is recommended.

The subsequent stages of application depend on the indication (see sections [A] - [G]).

When using light-curing materials, ORBI-Bond Universal can be applied without the addition of **ORBI-Bond Universal Dual Cure Activator** (see sections [A] - [E]).

The **ORBI-Bond Universal Dual Cure Activator** needs to be added only when working with self-curing or dual-curing materials (see sections [F] - [G]).

**[A] Direct or indirect restorations with light-curing restorative, core build-up or luting composites**

- Prepare the cavity and core according to the principles of the adhesive restoration technique, and remove any plaque or tartar in the case of noncarious class V cavities. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.

- Etching of the dental hard tissue (optional)
In self-etch mode, ORBI-Bond Universal achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is unprepared must always be etched at the beginning of the bonding process.
  - Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
  - Aspirate the acid-etch agent, rinse with water for approx. 20 s.
  - Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.
- Place one drop of ORBI-Bond Universal on a mixing palette.

- Apply the adhesive evenly to the surfaces of the cavity and rub it in for 20 s with a disposable applicator.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
- Apply the restorative, core build-up or luting material in accordance with the respective instructions for use.

**[B] Intraoral repair of fillings or restorations**

(made from composite, glass ceramic, zirconium oxide ceramic, aluminium oxide ceramic, metal (incl. amalgam))

- Thoroughly roughen the surface of the restoration to which the adhesive is to be applied and sandblast it, preferably with aluminium oxide < 50 µm. In the case of ceramics, remove the varnish 1 mm beyond the margin.
- Clean thoroughly with water spray and dry with moisture- and oil-free air.
- Continue with the application of ORBI-Bond Universal immediately afterwards (see [A] 3. - 6.). The use of an additional primer is not necessary.
- Carry out the repair using a light-curing composite according to the manufacturer’s instructions for use.

**[C] Treatment of hypersensitive tooth necks**

- Clean the tooth neck in the usual way (e.g. with polishing paste or using a powder jet).
- Clean thoroughly with water spray and dry with moisture- and oil-free air. Do not overdry. The objective is a surface with a silky matt appearance.
- Afterwards, continue with the application of ORBI-Bond Universal (see [A] 3. - 6.).
- The oxygen inhibition layer should be removed using a cloth or a foam pellet.

**[D] Protective varnish for glass ionomer cement restorations**

- If a matrix was used to facilitate shaping, remove the matrix first.
- Place one drop of ORBI-Bond Universal on a mixing palette.
- Apply ORBI-Bond Universal to the modelled filling with a disposable applicator.
- Dry off ORBI-Bond Universal with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure ORBI-Bond Universal on every surface of the restoration using a commercially available polymerisation device for 10 s (LED or light with an output of > 500 mW/cm²).
- Once the restoration has completely cured, grind the occlusion and seal the ground areas once more using ORBI-Bond Universal in order to prevent the restoration from drying out.
- The oxygen inhibition layer should be removed using a cloth or a foam pellet.

**[E] Sealing of cavities prior to amalgam restorations or prior to temporary luting**

Note: Do not seal cavities with ORBI-Bond Universal if temporary luting using a methacrylate-based material is scheduled afterwards.

- Prepare the cavity according to the further treatment planned. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
- Place one drop of ORBI-Bond Universal on a mixing palette.

- Apply ORBI-Bond Universal to the cavity surface to be treated with a disposable applicator and rub it in for 20 s.
- Dry off ORBI-Bond Universal with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure ORBI-Bond Universal using a commercially available polymerisation device for 10 s (LED or halogen light with an output of > 500 mW/cm²).
- Depending on the indication, continue to work with the desired material afterwards in accordance with the corresponding instructions for use.

**[F] Direct or indirect restorations with self- or dual-curing restorative, core build-up or luting composites**

In combination with ORBI-Bond Universal Dual Cure Activator, ORBI-Bond Universal creates a strong bond to dual-curing and self-curing methacrylate-based composites and luting materials.

- Prepare the cavity or core according to the principles of the adhesive restoration technique. Clean with water. Dry off excess water with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance. Protect areas near to the pulp with a suitable cavity lining based on calcium hydroxide.
- Etching of the dental hard tissue (optional)
In the self-etch mode, ORBI-Bond Universal achieves high levels of adhesive strength on enamel and dentine. Prior etching of the dental hard tissue may optionally be carried out. Only enamel which is unprepared must always be etched at the beginning of the bonding process (e.g. in the luting of Maryland bridges).
  - Etch the surfaces to be etched using a dental acid-etch agent based on phosphoric acid (see applicable instructions for use).
  - Aspirate the acid-etch agent, rinse with water for approx. 20 s.
  - Dry off excess moisture with a gentle stream of air. Do not overdry the dentine. The objective is a surface with a silky matt appearance.
- Mix one drop of ORBI-Bond Universal and one drop of ORBI-Bond Universal Dual Cure Activator in a mixing palette with a disposable applicator for 3 s.
- Apply the mixed adhesive evenly on the surfaces of the cavity and rub it in for 20 s with a disposable applicator.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents.
- Cure the adhesive layer for 10 s using a commercially available polymerisation device (LED or halogen light with an output of > 500 mW/cm²).
- Depending on the indication, continue to work with the desired material in accordance with the corresponding instructions for use.

**[G] Luting of root posts**

- Thoroughly mix one drop of ORBI-Bond Universal and one drop of ORBI-Bond Universal Dual Cure Activator in a mixing palette with a suitable disposable applicator (e.g. ORBI-Brush) for 3 s.
- Use the applicator to apply the mixed adhesive into the prepared root canal and to the occlusal surface of the core and rub it in for 20 s.
- Dry off the adhesive layer with dry, oil-free air for at least 5 s in order to remove any solvents. Remove any excess material by dabbing with an absorbent paper strip. Do not light-cure the adhesive layer.
- Introduce the root post, which has been prepared in accordance with the manufacturer’s instructions, into the root canal using a suitable self- or dualcuring composite for post luting (see applicable instructions for use) in such a way that a slight excess of material is achieved.
- To fix the root post and cure the occlusally-applied ORBI-Bond Universal, the material should be light-cured for 40 s.

**Warnings/precautionary measures:**

- Phenolic substances (e.g. temporary materials containing eugenol) impair polymerisation and therefore must not come into contact with ORBI-Bond Universal.
- ORBI-Bond Universal contains organic acids, (di)methacrylates, amines (DABE) and BHT. ORBI-Bond Universal should not be used in case of known hypersensitivities (allergies) to any of these ingredients.
- Avoid contact with the oral mucosa. Contact with ORBI-Bond Universal may cause whitening of the oral mucosa as a result of the coagulation of proteins. This is a reversible effect which subsides within a few days.
- In case of eye contact with ORBI-Bond Universal, rinse eyes immediately with plenty of water and seek medical advice.
- ORBI-Bond Universal is light-curing. Hence, please avoid exposure to strong lighting during application.
- The mixed adhesive is self-curing and light-curing, further processing must therefore take place immediately after mixing.

**Storage:**

Tightly close ORBI-Bond Universal and ORBI-Bond Universal Dual Cure Activator immediately after use. Store at temperatures of 2 °C to 8 °C (ORBI-Bond Universal) or rather 4 °C to 23 °C (ORBI-Bond Universal Dual Cure Activator).

Do not expose to direct sunlight or the operating light. Do not use after expiry date. Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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