ORBIS

ORBI-Sept SSC

Absauganlagenreiniger

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en

Method of administration:
Foam-free acidic cleaning concentrate for cleaning and care of dental suction systems. It is recommended to use the product twice per week.

Instructions for use:
The concentrate must first be diluted with water. Recommended application concentration: 5 %. For this purpose, top up 100ml of concentrate with tap water to 2000ml, suck it through the suction hose (saliva extractor) and leave it to take effect for at least 30 minutes. The additional cleaning via the mouthwash basin is carried out by pouring in the ready-to-use dilution.

Depending on the load, the extraction system should be cleaned once or twice a week with the suction system cleaner. Do not mix with other cleaners or disinfectants! Rinse the suction hoses and the cuspidor well with water before starting work. Only for professional users.

Safety instructions:
H315 Causes skin irritation. H 319 Causes severe eye irritation. Wear protective gloves/eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation occurs: Get medical advice/attention. If eye irritation persists: Get medical advice/attention. IF ON SKIN: Wash with plenty of water. Only for professional users.

(EC) 648/2004: < 5% amphoteric surfactants, < 5% phosphonates

All serious incidents related to the device shall be reported by the user/patient to the manufacturer and to the competent authority of the State in which the user/patient is established.

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